


botiss grafter

APPLICATOR FOR BONE SUBSTITUTE MATERIALS



SYSTEM THINKING:

ALLOWS APPLICATION OF VARIOUS BONE SUBSTITUTE MATERIALS OR THEIR MIXTURES

Autografts • Allografts • Bovine Bone • Synthetics

botiss grafter

Sieve cap enables optimal
hydration within the applicator

Efficient application of the desired amounts
of bone substitute materials

Application of up to 2 ml
of particulate bone substitute materials

Comfortable and clean handling

COMPACT **INTUITIVE DESIGN**



Provided sterile

Designed for **single use only**

Ideal for maxgraft[®], cerabone[®], maxresorb[®]


Allows **mixing of materials**




workflow




FILLING
via integrated funnel




INSERTION
of plunger rod




HYDRATION & EJECTION
of excess liquid
(avoid compression of granules)



REMOVAL of
SIEVE CAP



APPLICATION
of bone substitute material



Art.-No.	Content
550001	1 x botiss grafter



botiss biomaterials GmbH
Hauptstr. 28
15806 Zossen / Germany

Tel: +49 33769 / 88 41 985
Fax: +49 33769 / 88 41 986

www.botiss.com
[#botissdental](https://twitter.com/botissdental)



biotrics bioimplants AG, Ullsteinstr. 108, 12109 Berlin, Deutschland
www.biotrics.com

Rev.: BBGen-03/2024-01